

2 HANDS, 30 SECONDS KNOW MUCH MORE OF YOUR HEART



Handheld Heart Observer MD100E

Handheld Heart Observer MD100E helps users record the ECG waveform and know their heart condition via simple and quick measurement, so that they can improve their poor living habits, and better regulate their heart condition according to doctor's guidance.

CHOICEMMED™

Handheld Heart Observer MD100E

CVD is always threatening people's health

17,000,000 people died of heart attack annually!

3,000,000 people died as sudden death!

100,000,000 families are threatened by CVD!

A substantial way to survive from sudden death is the timely monitoring for arrhythmias! * CVD (cardiovascular disease)

2 HANDS, 30 SECONDS, KNOW MUCH MORE OF YOUR HEART

CHOICE Handheld Heart Observer, which is developed based on the heart disease "transient nature", simple and rapid measurement, small size, light weight and easy to carry.

CHOICE Handheld Heart Observer, to have the moment measurement in your family when you do not feel very well, to capture ECG waveforms timely when having heart disease incidence, to assist doctors to make early diagnosis.

We should change the passive status about our heart-health conditions, to have active control of our heart health.



A brilliant carry-on Heart Observer for family heart health



Enhance business man's confidence in success!



Enhance business man's confidence in success!



Life became more harmony!

Handheld Heart Observer MD100E

- * 2.8" large TFT screen is convenient senior citizen to read
- * 2G micro-SD card provides mass storage, information is easy to transfer
- * 250Hz sampling rate can capture subtle changes in ECG waveform
- * Low consumption: 2 AA alkaline batteries can be used not less than 500 measurements
- * Smiling face or desponding face will show your measurement result that is easy to know your heart statues
- * Analysis software will help doctor analyze ECG waveform subtle changes through adjusting display scale
- * 30 seconds and continuous measurement two modes: In the continuous measurement mode, HEART OBSERVER MD100E is a miniaturized single channel "electrocardiograph".



Brief Guide for Correct Measurement


How to use the device

Press the button to switch on the power, and then place your finger closely onto the finger electrodes.



Appropriate Posture

There are two options available. Posture One: Place the chest electrode against your left palm. Posture Two: Place the chest electrode against bare skin about 5 cm below the nipple.



Provide for better life

Beijing Choice Electronic Technology Co., Ltd

Address: Rm 1127-1128, Bailangyuan Building B, Fuxing Road A36, Beijing, P. R. China Zip: 100039
Tel: +86 10 88204188, 68221006, 88203520 FAX: +86 10 88204632
Email: market@choicemmed.com Website: www.choicemmed.com

